Identification Data



January 13, 2022

LAB GROWN DIAMOND Certificate No: 320060526

Gemprint is the unique optical fingerprint for positive identification of your lab grown diamond. Register your lab grown diamond at www.Gemprint.com and receive insurance discounts up to 10%.



Laser Inscription:

The illustration depicts enlarged and approximate appearances of the inscriptions. Girdle laser inscribed "GROWN IN THE USA BY WD", "PAT. 6,858,078", GCAL Logo and "LG320060526"







580 Fifth Ave LL-05 T 212-869-8985 GCALUSA.com



ISO/IEC 17025 2017 ANAB L2177-1 Accredited Testing Lab

The 4Cs Grading Analysis

GCAL 320060526 LAB GROWN DIAMOND*

Carat Weight: 1.56

Cut: Excellent **Oval Brilliant** Shape: Measurements: 8.98x6.64x4.13mm Optical Brilliance: Excellent Optical Symmetry: Very Good Polish: Excellent External Symmetry: Excellent Girdle Thickness: Medium-Thick Culet Size: None

G Color: Fluorescence: None

Clarity: Identifying Characteristic(s) Characteristic Location(s):

VS2 Crystals/Clouds Bezel, Table/Pavilion

*Comments: This laboratory grown diamond was created by the CVD (Chemical Vapor Deposition) method, and has the same chemical, physical, and optical properties as a mined diamond. This diamond is Type IIa, which means it is devoid of nitrogen impurities. As Grown - No evidence of post-growth treatment was detected.

Photomicrographs:

Actual images of the crown (top) and pavilion (bottom) of this diamond photographed at magnifications up to 10x.





Light Performance Profile

Optical Brilliance Analysis: Brilliance is the overall return of light to the viewer. The brilliance image is a representation of (a) white areas of light return, or brilliance, and (b) dark-blue areas of light loss.



Optical Brilliance

Optical Symmetry Analysis:

The colored areas of the symmetry image are indications of light handling ability, giving a visual representation of proportions and facet alignment.



Optical Symmetry

Proportion Diagram:

The proportion diagram illustrates the actual dimensions as recorded by optical scanning technology.

