Identification Data



October 8, 2021

LAB GROWN DIAMOND Certificate No: 312590340





Gemprint is the unique optical identification fingerprint of your lab grown diamond. Register your lab grown diamond fingerprint at www.Gemprint.com and receive insurance discounts up to 10%.


Laser Inscription



Girdle laser inscribed: GCAL LG312590340 GROWN IN THE USA BY WD PAT. 6,858,078 This illustration depicts the approximate appearance of the inscriptions



All certified certificate, ONLY available at an



GCALUSA.com



ANAB L2177-1 Accredited Testing Lab



The 4Cs Grading Analysis

GCAL 312590340 LAB GROWN DIAMOND*

Carat Weight: 0.76

Cut: Excellent **Oval Brilliant** Shape: Measurements: 7.00x5.18x3.28mm Optical Brilliance: Excellent Optical Symmetry: Very Good Polish: Very Good External Symmetry: Very Good Girdle Thickness: Medium-Thick Culet Size: None

Color: Fluorescence: None

Clarity: Identifying Characteristic(s) Characteristic Location(s):

Clouds/Crystal Upper Girdle, Table/Table

G

*Comments: This laboratory grown diamond was created by the CVD (Chemical Vapor Deposition) method, and has the same chemical, physical, and optical properties as a mined diamond. This diamond is Type IIa, which means it is devoid of nitrogen impurities. As Grown - No evidence of post-growth treatment was detected.

Photomicrographs:

Actual images of the crown (top) and pavilion (bottom) of this diamond photographed at magnifications up to 10x.

© 2021 GCAL





Light Performance Profile

Optical Brilliance Analysis: Brilliance is the overall return of light to the viewer. The brilliance image is a representation of (a) white areas of light return, or brilliance, and (b) dark-blue areas of light loss.



Optical Brilliance

Optical Symmetry Analysis:

The colored areas of the symmetry image are indications of light handling ability, giving a visual representation of proportions and facet alignment.



Optical Symmetry

Proportion Diagram:

The proportion diagram illustrates the actual dimensions as recorded by optical scanning technology.

